

People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively

FM made friendly
- and cost-effective

child friendly hearing care

It takes a truly dedicated approach to help children with hearing problems achieve their full potential. That's why we deliver the solutions and services that professionals and caregivers need to provide children the opportunities they deserve. This is what child-friendly hearing care is all about.



124855UK/05.12 Printed on 100% recycled paper



child friendly hearing care



The **flexible** and cost **effective** FM solution

With the Amigo range of FM equipment, Oticon offers a quality FM solution in all price categories.

The Amigo family is characterised by its reliability and ease of use, and the T5 carries on this unique tradition.

Certainty The LED light on the Amigo T5 provides instant visual confirmation that the transmitter is fully charged, charging, or muted.

Durability The internal components and external casing have been designed to withstand continuous use in demanding environments such as schools.

Flexibility The Amigo T5 is fully compatible with most FM receivers and transmitters, cochlear implants and hearing aids.

Convenience The operation and set-up of the Amigo T5 is extremely user friendly. The easy-to-use push buttons combined with a clear display make system use and troubleshooting very easy.

Cost effective The T5 represents the most cost efficient solution in the Amigo family. So even more children can gain access to clear sound!

Amigo T5 transmitter

- External antenna built into the microphone cord
- One-touch mute
- Fail safe charging
- FM transmission range ~ 30 m
- LED light



How **Amigo** can help

Fail-safe battery charging

The Amigo T5 uses rechargeable AAA batteries, but it is safe to use standard alkaline batteries because of the intelligent battery detection scheme.

Quick & Easy Mute


The microphone can be muted when the FM signal is of no actual value - e.g. when the teacher is moving around in the classroom but not actually speaking.

Clear and simple user controls

The informative display, clearly marked jacks and easy to operate buttons make the Amigo T5 user friendly, even for inexperienced users.

LED on the Amigo

The built-in LED light let you know that the system is charging, the battery level is low or the transmitter is muted.



Amigo T5 has three different microphone options depending on your need:
- Directional lapel
- Omni directional lapel
- Boom



Amigo T5 comes with a wall charger that can charge 2 units simultaneously.



The Amigo T5 can be connected to external audio sources (such as a TV, PC or MP3 player) with the aid of auxiliary cables.

